



12 Lead Portable ECG

DECK MOUNT[®]
ELECTRONICS

Rhythm Ultra Max

Now capture **ECG** in just 10 seconds



Portable ECG with 24-bit ADC precision, offering full cardiac signal capture, secure cloud connectivity, instant reporting, and clinical-grade accuracy


ADVANCED FEATURES



FLUKE
TESTED

Engineered for Clinical Excellence

Rhythm Ultra Max combines cutting-edge technology with practical design to deliver accurate, reliable, and accessible cardiac diagnostics anywhere, anytime.



Hospital-Grade 12-Lead Accuracy

Simultaneous 12-lead ECG from 10 electrodes, tested on **Fluke gold-standard** equipment, with 1000 samples/sec/channel for precise and reliable cardiac assessment.



Extended Battery Performance

4000 mAh battery delivers up to 1 month of standby or 4hrs backup on a single charge. Reliable performance with minimal downtime and freedom from frequent recharging.



Portable Handheld Design

Weighing just 298 grams, this lightweight, handheld ECG is easy to carry and ideal for on-the-go, bedside, and remote care.



Rhythm Ultra Max

Mobile App-Based ECG Recording

Record ECGs seamlessly through a mobile app with instant data capture, real-time viewing, and effortless report sharing.



99.9%
ACCURACY

99.9% Accuracy

Delivering exceptional precision for confident diagnosis and improved patient outcomes.



CardioX

Smart ECG Monitoring Ecosystem

An ECG monitoring system delivering live cardiac analytics and clinical insights built for India's healthcare ecosystem.



Live ECG Monitoring

Real-time cardiac waveform capture and display with clinical-grade 12-lead ECG output across hospital, clinic, and remote settings.



Hyperkalemia + HRV Detection

Flags Hyperkalemia and Hypokalemia via T-wave analysis. HRV metrics (SDNN, RMSSD) assessed automatically on every recording.



Analytics Dashboard

A unified digital dashboard for clinicians view trends, patient history, and cardiac reports in one integrated interface.



1. Add Doctors



2. Share ECG Reports



3. Track Review Status

Date	Time	Op	Doctor	Test No	Age	Gender	Height	Weight	Type	Status
2024-05-07	15:40:46								ECG	Pending
2024-05-07	15:40:46								ECG	Pending
2024-05-07	15:40:55								ECG	Pending
2024-05-07	15:40:55								ECG	Pending
2024-05-07	15:41:11								ECG	Pending
2024-05-07	15:41:11								ECG	Pending
2024-05-07	15:41:22								ECG	Pending
2024-05-07	15:41:22								ECG	Pending
2024-05-07	15:41:36								ECG	Pending
2024-05-07	15:41:36								ECG	Pending
2024-05-07	15:41:36								ECG	Pending
2024-05-08	11:38:53								ECG	Pending
2024-05-08	11:38:53								ECG	Pending
2024-05-08	11:38:53								ECG	Pending
2024-05-08	11:38:53								ECG	Pending
2024-05-08	11:38:53								ECG	Pending



Complete 12-Lead Dashboard

All 12 leads visible simultaneously with axis analysis, waveform overlay, and real-time signal quality indicator.



Multi-Format ECG Reports

Generate reports in 4x3, 12x1, or 6x2 layouts in US and Cabrera (European) standards – share or print instantly.



Doctor Database & Sharing

Maintain a doctor directory, send reports directly, track review status, and receive "Reviewed by" confirmation on shared ECGs.



Hospitals & ICUs



Clinics & OPDs



Remote & Home Care

Product Specifications

Product Name	RHYTHM ULTRA MAX
Dimensions	171 mm (L) × 95 mm (W) × 40.5 mm (H)
Weight	298 grams
Power	Rechargeable Li-Po Battery, 3.7 V/ 4000 mAh
Charging	USB Type B, 6-8 hrs Charging Time 1 Month Standby or Up to 160 ECGs
Operating Conditions	10 °C to 50 °C (50 °F to 122°F) 20% to 90% Relative Humidity
Storage Conditions	0 °C to 55 °C (32 °F to 131°F) 10% to 90% Relative Humidity
Warranty	1 year

Patient Data	Name, Age, Gender
LED Indicators	LED indication on the device for the battery, Electrodes, App Connection
Analysis	- RR, HR, PR, QRS, QT, QTc, QTcf, RV5/SV1, RV5+SV1 P/QRS/T (Deviation)
ECG Sharing (PDF)	WhatsApp™, Mail
App Compatibility	Supports Windows™10 & above, Android™9.0 & above

What's in the Box	Rhythm Ultra Max, Power cable and adapter, Electrodes x10, Cable with Crocodile clip/ Banana cable with clamps (optional), Carry bag, User manual
-------------------	--

ECG Acquisition Details

Number of Channels	12-lead Simultaneous ECG from 10 electrodes
Connectivity	BLE 5.0(Bluetooth™) & USB 2.0 & above
Frequency Response	0.01 ~ 250 Hz
Time Constant	> 3 Seconds
Input Impedance	> 50M Ohms
Input CIR Current	≤ 0.1 μA
Patient Leak Current	< 10 μA
Sensitivity (Amplitude)	5, 10, 20 mm/ mV
Recording Speed	12.5, 25, 50 mm/ Sec
CMRR	> 110 dB
A to D Conversion	24 bit
Sampling Frequency	1000 Samples per sec/ channel
Filters	Line Frequency, Base Line, and Muscle Filters
Noise Level	< 15μVp-p



Product Code : PEEG

Manufactured by

Deck Mount Electronics Private Limited

683, Udyog Vihar phase V, Sector 19, Gurugram. Haryana 122016

